

产品名称: CD49B Mouse Monoclonal Antibody
产品货号: AMM82060



产品概述 (Summary)

产品名称 (Production Name)	CD49B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD49B
别名 (Alternative Names)	ITGA2; BR; GPIa; HPA-5; VLA-2; VLAA2
基因 ID (Gene ID)	3673.0
蛋白 ID (SwissProt ID)	P17301.Purified recombinant fragment of human CD49B (AA: extra 994-1132) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:100-1:500,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	129.3kDa

研究背景 (Background)

This gene encodes the alpha subunit of a transmembrane receptor for collagens and related proteins. The encoded protein

产品名称: CD49B Mouse Monoclonal Antibody
产品货号: AMM82060




forms a heterodimer with a beta subunit and mediates the adhesion of platelets and other cell types to the extracellular matrix. Loss of the encoded protein is associated with bleeding disorder platelet-type 9. Antibodies against this protein are found in several immune disorders, including neonatal alloimmune thrombocytopenia. This gene is located adjacent to a related alpha subunit gene. Alternative splicing results in multiple transcript variants.


研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD49B mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: CD49B Mouse Monoclonal Antibody
产品货号: AMM82060



Flow cytometric analysis of HT1080 cells using CD49B mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .