

产品名称: CD335 Mouse Monoclonal Antibody
产品货号: AMM82068



产品概述 (Summary)

产品名称 (Production Name)	CD335 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD335
别名 (Alternative Names)	NCR1; LY94; NKP46; NK-p46
基因 ID (Gene ID)	9437.0
蛋白 ID (SwissProt ID)	O76036.Purified recombinant fragment of human CD335 (AA: extra 22-258) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	34.5kDa

研究背景 (Background)

NCR1 (Natural Cytotoxicity Triggering Receptor 1) is a Protein Coding gene. Diseases associated with NCR1 include Viral


产品名称: CD335 Mouse Monoclonal Antibody
产品货号: AMM82068




Infectious Disease. Among its related pathways are Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell and Natural Killer Cell Receptors: Human Target Cell – NK Cell Ligand-Receptor Interactions. GO annotations related to this gene include signal transducer activity, downstream of receptor. An important paralog of this gene is LILRB3.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD335 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .