

产品名称: CEACAM8 Mouse Monoclonal Antibody
产品货号: AMM82070



产品概述 (Summary)

产品名称 (Production Name)	CEACAM8 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEACAM8
别名 (Alternative Names)	CEACAM8; CD67; CGM6; NCA-95
基因 ID (Gene ID)	1088.0
蛋白 ID (SwissProt ID)	P31997.Purified recombinant fragment of human CEACAM8 (AA: 43-178) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	38.2kDa

研究背景 (Background)

CEACAM8 (Carcinoembryonic Antigen Related Cell Adhesion Molecule 8) is a Protein Coding gene. Diseases associated


产品名称: CEACAM8 Mouse Monoclonal Antibody
产品货号: AMM82070




with CEACAM8 include Endocarditis and Shwartzman Phenomenon. Among its related pathways are Innate Immune System and Cell surface interactions at the vascular wall. An important paralog of this gene is CEACAM1.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CEACAM8 mouse mAb against Hela (1) cell lysate.

注意事项 (Note)

For research use only .