

产品名称: CD64 Mouse Monoclonal Antibody
产品货号: AMM82073



产品概述 (Summary)

产品名称 (Production Name)	CD64 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD64
别名 (Alternative Names)	FCGR1A; FCRI; CD64A; IGFR1
基因 ID (Gene ID)	2209.0
蛋白 ID (SwissProt ID)	P12314.Purified recombinant fragment of human CD64 (AA: extra 16-145) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	42.6kDa

研究背景 (Background)

This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma


产品名称: CD64 Mouse Monoclonal Antibody
产品货号: AMM82073




receptor. The gene is one of three related gene family members located on chromosome 1.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD64 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using CD64 mouse mAb with DAB staining.

产品名称: CD64 Mouse Monoclonal Antibody
产品货号: AMM82073



注意事项 (Note)

For research use only .