

产品名称: CHRNE Mouse Monoclonal Antibody
产品货号: AMM82075



产品概述 (Summary)

产品名称 (Production Name)	CHRNE Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHRNE
别名 (Alternative Names)	ACHRE; CMS1D; CMS1E; CMS2A; CMS4A; CMS4B; CMS4C; FCCMS; SCCMS
基因 ID (Gene ID)	1145.0
蛋白 ID (SwissProt ID)	Q04844.Purified recombinant fragment of human CHRNE (AA: extra 21-239) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	54.7kDa

研究背景 (Background)

Acetylcholine receptors at mature mammalian neuromuscular junctions are pentameric protein complexes composed of


产品名称: CHRNE Mouse Monoclonal Antibody
产品货号: AMM82075




four subunits in the ratio of two alpha subunits to one beta, one epsilon, and one delta subunit. The acetylcholine receptor changes subunit composition shortly after birth when the epsilon subunit replaces the gamma subunit seen in embryonic receptors. Mutations in the epsilon subunit are associated with congenital myasthenic syndrome.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CHRNE mouse mAb against C6 (1) cell lysate.



Flow cytometric analysis of SK-N-SH cells using CHRNE mouse mAb (green) and negative control (red).

产品名称: CHRNE Mouse Monoclonal Antibody
产品货号: AMM82075



注意事项 (Note)

For research use only .