

产品名称: TNFRSF10D Mouse Monoclonal Antibody
产品货号: AMM82089



产品概述 (Summary)

产品名称 (Production Name)	TNFRSF10D Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFRSF10D
别名 (Alternative Names)	DCR2; CD264; TRUNDD; TRAILR4; TRAIL-R4
基因 ID (Gene ID)	8793.0
蛋白 ID (SwissProt ID)	Q9UBN6.Purified recombinant fragment of human TNFRSF10D (AA: extra 56-211) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41.8kDa

研究背景 (Background)

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular

产品名称: TNFRSF10D Mouse Monoclonal Antibody
产品货号: AMM82089




TRAIL-binding domain, a transmembrane domain, and a truncated cytoplasmic death domain. This receptor does not induce apoptosis, and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis.


研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using TNFRSF10D mouse mAb against Hela (1), L-02 (2), A549 (3), HepG2 (4), and COS7 (5) cell lysate.

注意事项 (Note)

For research use only .