

产品名称: CD201 Mouse Monoclonal Antibody
产品货号: AMM82093



产品概述 (Summary)

产品名称 (Production Name)	CD201 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD201
别名 (Alternative Names)	CCCA; EPCR; CCD41
基因 ID (Gene ID)	10544.0
蛋白 ID (SwissProt ID)	Q9UNN8.Purified recombinant fragment of human CD201 (AA: extra 18-210) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	26.7kDa

研究背景 (Background)

The protein encoded by this gene is a receptor for activated protein C, a serine protease activated by and involved in the


产品名称: CD201 Mouse Monoclonal Antibody
产品货号: AMM82093



blood coagulation pathway. The encoded protein is an N-glycosylated type I membrane protein that enhances the activation of protein C. Mutations in this gene have been associated with venous thromboembolism and myocardial infarction, as well as with late fetal loss during pregnancy. The encoded protein may also play a role in malarial infection and has been associated with cancer.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)

注意事项 (Note)

For research use only .