

产品名称: LILRA3 Mouse Monoclonal Antibody
产品货号: AMM82101



产品概述 (Summary)

产品名称 (Production Name)	LILRA3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRA3
别名 (Alternative Names)	HM31; HM43; ILT6; LIR4; CD85E; ILT-6; LIR-4
基因 ID (Gene ID)	11026.0
蛋白 ID (SwissProt ID)	Q8N6C8.Purified recombinant fragment of human LILRA3 (AA: 24-168) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	47.5kDa

研究背景 (Background)

This gene encodes a member of a family of immunoreceptors that are expressed predominantly in monocytes and B cells,


产品名称: LILRA3 Mouse Monoclonal Antibody
产品货号: AMM82101




and at lower levels in dendritic cells and natural killer cells. The encoded protein lacks the transmembrane region found in other members of this family. It acts as a soluble receptor for class I major histocompatibility complex (MHC) antigens. Alternatively spliced transcript variants encoding different isoforms have been found. This gene is located in a cluster of related genes on chromosome 19 and is polymorphic in human populations, with many individuals containing a deletion of this genomic region.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of K562 cells using LILRA3 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .