

产品名称: CD46 Mouse Monoclonal Antibody
产品货号: AMM82106



产品概述 (Summary)

产品名称 (Production Name)	CD46 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD46
别名 (Alternative Names)	MCP; TLX; AHUS2; MIC10; TRA2.10
基因 ID (Gene ID)	4179.0
蛋白 ID (SwissProt ID)	P15529.Purified recombinant fragment of human CD46 (AA: extra 35-179) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ICC 1:100-1:500,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	43.7kDa

研究背景 (Background)

The protein encoded by this gene is a type I membrane protein and is a regulatory part of the complement system. The


产品名称: CD46 Mouse Monoclonal Antibody
产品货号: AMM82106




encoded protein has cofactor activity for inactivation of complement components C3b and C4b by serum factor I, which protects the host cell from damage by complement. In addition, the encoded protein can act as a receptor for the Edmonston strain of measles virus, human herpesvirus-6, and type IV pili of pathogenic *Neisseria*. Finally, the protein encoded by this gene may be involved in the fusion of the spermatozoa with the oocyte during fertilization. Mutations at this locus have been associated with susceptibility to hemolytic uremic syndrome. Alternatively spliced transcript variants encoding different isoforms have been described.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CD46 mouse mAb against SW480 (1) and MCF-7 (2) cell lysate.

产品名称: CD46 Mouse Monoclonal Antibody
产品货号: AMM82106



Immunofluorescence analysis of HeLa cells using CD46 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HL-60 cells using CD46 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .