

产品名称: TNFRSF19 Mouse Monoclonal Antibody
产品货号: AMM82109



产品概述 (Summary)

产品名称 (Production Name)	TNFRSF19 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFRSF19
别名 (Alternative Names)	TAJ; TROY; TRADE; TAJ-alpha
基因 ID (Gene ID)	55504.0
蛋白 ID (SwissProt ID)	Q9NS68 .Purified recombinant fragment of human TNFRSF19 (AA: extra 30-170) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	46kDa

研究背景 (Background)

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is highly expressed during


产品名称: TNFRSF19 Mouse Monoclonal Antibody
产品货号: AMM82109




embryonic development. It has been shown to interact with TRAF family members, and to activate JNK signaling pathway when overexpressed in cells. This receptor is capable of inducing apoptosis by a caspase-independent mechanism, and it is thought to play an essential role in embryonic development. Alternatively spliced transcript variants encoding distinct isoforms have been described.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using TNFRSF19 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .