

产品名称: CD360 Mouse Monoclonal Antibody  
产品货号: AMM82125



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | CD360 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody       |
| 宿主 (Host)              | Mouse                           |
| 应用 (Application)       | ELISA,FC                        |
| 种属反应性 (Reactivity)     | Human                           |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                              |
|------------------------|----------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CD360                                                                                        |
| 别名 (Alternative Names) | IL21R; NILR; IMD56                                                                           |
| 基因 ID (Gene ID)        | 50615.0                                                                                      |
| 蛋白 ID (SwissProt ID)   | Q9HBE5.Purified recombinant fragment of human CD360 (AA: extra 20-232) expressed in E. Coli. |

## 产品应用 (Application)

|                          |                                     |
|--------------------------|-------------------------------------|
| 稀释比 (Dilution Ratio)     | ELISA 1:5000-1:20000,FC 1:200-1:400 |
| 蛋白分子量 (Molecular Weight) | 59.1kDa                             |

## 研究背景 (Background)

The protein encoded by this gene is a cytokine receptor for interleukin 21 (IL21). It belongs to the type I cytokine receptors,

产品名称: CD360 Mouse Monoclonal Antibody  
产品货号: AMM82125




and has been shown to form a heterodimeric receptor complex with the common gamma-chain, a receptor subunit also shared by the receptors for interleukin 2, 4, 7, 9, and 15. This receptor transduces the growth promoting signal of IL21, and is important for the proliferation and differentiation of T cells, B cells, and natural killer (NK) cells. The ligand binding of this receptor leads to the activation of multiple downstream signaling molecules, including JAK1, JAK3, STAT1, and STAT3. Knockout studies of a similar gene in mouse suggest a role for this gene in regulating immunoglobulin production. Three alternatively spliced transcript variants have been described.


## 研究领域 (Research Area)

Jak-STAT signaling pathway

## 图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of K562 cells using CD360 mouse mAb (green) and negative control (red).

## 注意事项 (Note)

For research use only .