

产品名称: CD53 Mouse Monoclonal Antibody
产品货号: AMM82127



产品概述 (Summary)

产品名称 (Production Name)	CD53 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD53
别名 (Alternative Names)	MOX44; TSPAN25
基因 ID (Gene ID)	963.0
蛋白 ID (SwissProt ID)	P19397.Purified recombinant fragment of human CD53 (AA: extra mix) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	24.3kDa

研究背景 (Background)

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family.


产品名称: CD53 Mouse Monoclonal Antibody
产品货号: AMM82127




Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins. It contributes to the transduction of CD2-generated signals in T cells and natural killer cells and has been suggested to play a role in growth regulation. Familial deficiency of this gene has been linked to an immunodeficiency associated with recurrent infectious diseases caused by bacteria, fungi and viruses. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using CD53 mouse mAb against Raji (1), Ramos (2), Jurkat (3), MOLT4 (4), K562 (5), HL-60 (6), THP-1 (7), and U937 (8) cell lysate.


产品名称: CD53 Mouse Monoclonal Antibody
产品货号: AMM82127



Flow cytometric analysis of HL-60 cells using CD53 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD53 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using CD53 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .