

产品名称: CD353 Mouse Monoclonal Antibody
产品货号: AMM82140



产品概述 (Summary)

产品名称 (Production Name)	CD353 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD353
别名 (Alternative Names)	SLAMF8; BLAME; SBBI42
基因 ID (Gene ID)	56833.0
蛋白 ID (SwissProt ID)	Q9P0V8.Purified recombinant fragment of human CD353 (AA: extra 23-233) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	31.7kDa

研究背景 (Background)

This gene encodes a member of the CD2 family of cell surface proteins involved in lymphocyte activation. These proteins


产品名称: CD353 Mouse Monoclonal Antibody
产品货号: AMM82140




are characterized by Ig domains. This protein is expressed in lymphoid tissues, and studies of a similar protein in mouse suggest that it may function during B cell lineage commitment. The gene is found in a region of chromosome 1 containing many CD2 genes.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD353 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .