

产品名称: LILRB5 Mouse Monoclonal Antibody
产品货号: AMM82155



产品概述 (Summary)

产品名称 (Production Name)	LILRB5 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRB5
别名 (Alternative Names)	LIR8; CD85C; LIR-8
基因 ID (Gene ID)	10990.0
蛋白 ID (SwissProt ID)	O75023 .Purified recombinant fragment of human LILRB5 (AA: extra 24-164) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	64.1kDa

研究背景 (Background)

This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at


产品名称: LILRB5 Mouse Monoclonal Antibody
产品货号: AMM82155




chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). Several other LIR subfamily B receptors are expressed on immune cells where they bind to MHC class I molecules on antigen-presenting cells and inhibit stimulation of an immune response. Multiple transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using LILRB5 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .