

产品名称: RAB27B Mouse Monoclonal Antibody
产品货号: AMM82160



产品概述 (Summary)

产品名称 (Production Name)	RAB27B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB27B
别名 (Alternative Names)	C25KG
基因 ID (Gene ID)	5874.0
蛋白 ID (SwissProt ID)	O00194.Purified recombinant fragment of human RAB27B (AA: 2-218) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	24.6kDa

研究背景 (Background)

Members of the Rab protein family, including RAB27B, are prenylated, membrane-bound proteins involved in vesicular


产品名称: RAB27B Mouse Monoclonal Antibody
产品货号: AMM82160




fusion and trafficking (Chen et al., 1997 [PubMed 9066979]).

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using RAB27B mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebellum tissues using RAB27B mouse mAb with DAB staining.

产品名称: RAB27B Mouse Monoclonal Antibody
产品货号: AMM82160



注意事项 (Note)

For research use only .