

产品名称: CD141 Mouse Monoclonal Antibody
产品货号: AMM82163



产品概述 (Summary)

产品名称 (Production Name)	CD141 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD141
别名 (Alternative Names)	THBD; TM; THRM; AHUS6; BDCA3; THPH12
基因 ID (Gene ID)	7056.0
蛋白 ID (SwissProt ID)	P07204.Purified recombinant fragment of human CD141 (AA: extra 297-505) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	60.3kDa

研究背景 (Background)

The protein encoded by this intronless gene is an endothelial-specific type I membrane receptor that binds thrombin. This


产品名称: CD141 Mouse Monoclonal Antibody
产品货号: AMM82163




binding results in the activation of protein C, which degrades clotting factors Va and VIIIa and reduces the amount of thrombin generated. Mutations in this gene are a cause of thromboembolic disease, also known as inherited thrombophilia.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD141 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .