

产品名称: LRP6 Mouse Monoclonal Antibody
产品货号: AMM82170



产品概述 (Summary)

产品名称 (Production Name)	LRP6 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRP6
别名 (Alternative Names)	ADCAD2; STHAG7
基因 ID (Gene ID)	4040.0
蛋白 ID (SwissProt ID)	O75581.Purified recombinant fragment of human LRP6 (AA: extra 1198-1370) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	180kDa

研究背景 (Background)

This gene encodes a member of the low density lipoprotein (LDL) receptor gene family. LDL receptors are transmembrane

产品名称: LRP6 Mouse Monoclonal Antibody
产品货号: AMM82170




cell surface proteins involved in receptor-mediated endocytosis of lipoprotein and protein ligands. The protein encoded by this gene functions as a receptor or, with Frizzled, a co-receptor for Wnt and thereby transmits the canonical Wnt/beta-catenin signaling cascade. Through its interaction with the Wnt/beta-catenin signaling cascade this gene plays a role in the regulation of cell differentiation, proliferation, and migration and the development of many cancer types. This protein undergoes gamma-secretase dependent RIP- (regulated intramembrane proteolysis) processing but the precise locations of the cleavage sites have not been determined.


研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using LRP6 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .