

产品名称: CD49E Mouse Monoclonal Antibody
产品货号: AMM82195



产品概述 (Summary)

产品名称 (Production Name)	CD49E Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD49E
别名 (Alternative Names)	ITGA5; FNRA; VLA-5; VLA5A
基因 ID (Gene ID)	3678.0
蛋白 ID (SwissProt ID)	P08648.Purified recombinant fragment of human CD49E (AA: extra 111-253) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	114.5kDa

研究背景 (Background)

The product of this gene belongs to the integrin alpha chain family. Integrins are heterodimeric integral membrane

产品名称: CD49E Mouse Monoclonal Antibody
产品货号: AMM82195




proteins composed of an alpha subunit and a beta subunit that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha 5 subunit. This subunit associates with the beta 1 subunit to form a fibronectin receptor. This integrin may promote tumor invasion, and higher expression of this gene may be correlated with shorter survival time in lung cancer patients. Note that the integrin alpha 5 and integrin alpha V subunits are encoded by distinct genes.


研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD49E mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .