

产品名称: CD192 Mouse Monoclonal Antibody
产品货号: AMM82211



产品概述 (Summary)

产品名称 (Production Name)	CD192 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD192
别名 (Alternative Names)	CCR2; CKR2; CCR-2; CCR2A; CCR2B; CKR2A; CKR2B; CMKBR2; MCP-1-R; CC-CKR-2
基因 ID (Gene ID)	729230.0
蛋白 ID (SwissProt ID)	P41597 .Purified recombinant fragment of human CD192 expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	41.9kDa

研究背景 (Background)

The protein encoded by this gene is a receptor for monocyte chemoattractant protein-1, a chemokine which specifically


产品名称: CD192 Mouse Monoclonal Antibody
产品货号: AMM82211




mediates monocyte chemotaxis. Monocyte chemoattractant protein-1 is involved in monocyte infiltration in inflammatory diseases such as rheumatoid arthritis as well as in the inflammatory response against tumors. The encoded protein mediates agonist-dependent calcium mobilization and inhibition of adenylyl cyclase. This protein can also be a coreceptor with CD4 for HIV-1 infection. This gene is located in the chemokine receptor gene cluster region of chromosome 3.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD192 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .