

产品名称: CD275 Mouse Monoclonal Antibody
产品货号: AMM82214



产品概述 (Summary)

产品名称 (Production Name)	CD275 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD275
别名 (Alternative Names)	ICOSLG; B7H2; GL50; B7-H2; B7RP1; ICOSL; LICOS; B7RP-1; ICOS-L
基因 ID (Gene ID)	23308.0
蛋白 ID (SwissProt ID)	O75144.Purified recombinant fragment of human CD275 (AA: extra 19-256) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	33.3kDa

研究背景 (Background)

Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-cell proliferation and cytokine


产品名称: CD275 Mouse Monoclonal Antibody
产品货号: AMM82214




secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-stimulating memory T-cell function (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD275 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .