

产品名称: CD27 Mouse Monoclonal Antibody
产品货号: AMM82219



产品概述 (Summary)

产品名称 (Production Name)	CD27 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD27
别名 (Alternative Names)	T14; S152; Tp55; TNFRSF7; S152. LPFS2
基因 ID (Gene ID)	939.0
蛋白 ID (SwissProt ID)	P26842.Purified recombinant fragment of human CD27 (AA: extra 20-144) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	29.1kDa

研究背景 (Background)

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is required for generation

产品名称: CD27 Mouse Monoclonal Antibody
产品货号: AMM82219




and long-term maintenance of T cell immunity. It binds to ligand CD70, and plays a key role in regulating B-cell activation and immunoglobulin synthesis. This receptor transduces signals that lead to the activation of NF-kappaB and MAPK8/JNK. Adaptor proteins TRAF2 and TRAF5 have been shown to mediate the signaling process of this receptor. CD27-binding protein (SIVA), a proapoptotic protein, can bind to this receptor and is thought to play an important role in the apoptosis induced by this receptor.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD27 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .