

产品名称: LILRA1 Mouse Monoclonal Antibody
产品货号: AMM82230



产品概述 (Summary)

产品名称 (Production Name)	LILRA1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRA1
别名 (Alternative Names)	LIR6; CD85I; LIR-6
基因 ID (Gene ID)	11024.0
蛋白 ID (SwissProt ID)	O75019 .Purified recombinant fragment of human LILRA1 (AA: extra 322-461) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	53.3kDa

研究背景 (Background)

This gene encodes an activating member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a


产品名称: LILRA1 Mouse Monoclonal Antibody
产品货号: AMM82230




gene cluster at chromosomal region 19q13.4. The encoded protein is predominantly expressed in B cells, interacts with major histocompatibility complex class I ligands, and contributes to the regulation of immune responses. Alternative splicing results in multiple transcript variants encoding different isoforms.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of MOLT4 cells using LILRA1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .