

产品名称: LRP4 Mouse Monoclonal Antibody
产品货号: AMM82234



产品概述 (Summary)

产品名称 (Production Name)	LRP4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LRP4
别名 (Alternative Names)	CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2
基因 ID (Gene ID)	4038.0
蛋白 ID (SwissProt ID)	O75096.Purified recombinant fragment of human LRP4 (AA: extra 21-189) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	212kDa

研究背景 (Background)

This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a

产品名称: LRP4 Mouse Monoclonal Antibody
产品货号: AMM82234




regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome. [provided by RefSeq, May 2010]


研究领域 (Research Area)

Wnt signaling pathway

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using LRP4 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .