

产品名称: CD307E Mouse Monoclonal Antibody
产品货号: AMM82247



产品概述 (Summary)

产品名称 (Production Name)	CD307E Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ICC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD307E
别名 (Alternative Names)	FCRL5; CD307; FCRH5; IRTA2; BXMAS1; PRO820
基因 ID (Gene ID)	83416.0
蛋白 ID (SwissProt ID)	Q96RD9.Purified recombinant fragment of human CD307E (AA: extra 16-158) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	106.4kDa

研究背景 (Background)

This gene encodes a member of the immunoglobulin receptor superfamily and the Fc-receptor like family. This gene and


产品名称: CD307E Mouse Monoclonal Antibody
产品货号: AMM82247




several other Fc receptor-like gene members are clustered on the long arm of chromosome 1. The encoded protein is a single-pass type I membrane protein and contains 8 immunoglobulin-like C2-type domains. This gene is implicated in B cell development and lymphomagenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD307E mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

产品名称: CD307E Mouse Monoclonal Antibody
产品货号: AMM82247



Flow cytometric analysis of HL-60 cells using CD307E mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .