

产品名称: MTA1 Mouse Monoclonal Antibody
产品货号: AMM82255



产品概述 (Summary)

产品名称 (Production Name)	MTA1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MTA1
别名 (Alternative Names)	MTA1
基因 ID (Gene ID)	9112.0
蛋白 ID (SwissProt ID)	Q13330.Purified recombinant fragment of human MTA1 (AA: 541-715) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	80.8kDa

研究背景 (Background)

This gene encodes a protein that was identified in a screen for genes expressed in metastatic cells, specifically, mammary


产品名称: MTA1 Mouse Monoclonal Antibody
产品货号: AMM82255




adenocarcinoma cell lines. Expression of this gene has been correlated with the metastatic potential of at least two types of carcinomas although it is also expressed in many normal tissues. The role it plays in metastasis is unclear. It was initially thought to be the 70kD component of a nucleosome remodeling deacetylase complex, NuRD, but it is more likely that this component is a different but very similar protein. These two proteins are so closely related, though, that they share the same types of domains. These domains include two DNA binding domains, a dimerization domain, and a domain commonly found in proteins that methylate DNA. The profile and activity of this gene product suggest that it is involved in regulating transcription and that this may be accomplished by chromatin remodeling. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2011]

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using MTA1 mouse mAb against SW480 (1), T47D (2), MCF-7 (3), and COS7 (4) cell lysate.

产品名称: MTA1 Mouse Monoclonal Antibody
产品货号: AMM82255



Flow cytometric analysis of Hela cells using MTA1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .