

产品名称: CD218A Mouse Monoclonal Antibody
产品货号: AMM82275



产品概述 (Summary)

产品名称 (Production Name)	CD218A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD218A
别名 (Alternative Names)	IL18R1; IL18RA; IL1RRP; CDw218a; IL-1Rrp
基因 ID (Gene ID)	8809.0
蛋白 ID (SwissProt ID)	Q13478 .Purified recombinant fragment of human CD218A (AA: extra 22-185) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	62.3kDa

研究背景 (Background)

The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor


产品名称: CD218A Mouse Monoclonal Antibody
产品货号: AMM82275




specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD218A mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .