

产品名称: RIPK1 Mouse Monoclonal Antibody
产品货号: AMM82284



产品概述 (Summary)

产品名称 (Production Name)	RIPK1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIPK1
别名 (Alternative Names)	RIP; RIP1; IMD57; RIP-1
基因 ID (Gene ID)	8737.0
蛋白 ID (SwissProt ID)	Q13546.Purified recombinant fragment of human RIPK1 (AA: 131-277) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	76kDa

研究背景 (Background)

This gene encodes a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases. The

产品名称: RIPK1 Mouse Monoclonal Antibody
产品货号: AMM82284




encoded protein plays a role in inflammation and cell death in response to tissue damage, pathogen recognition, and as part of developmental regulation. RIPK1/RIPK3 kinase-mediated necrosis is referred to as necroptosis. Genetic disruption of this gene in mice results in death shortly after birth. [provided by RefSeq, Aug 2017]


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using RIPK1 mouse mAb against HL-60 (1), K562 (2), and Jurkat (3) cell lysate.

注意事项 (Note)

For research use only .