

产品名称: CD315 Mouse Monoclonal Antibody
产品货号: AMM82296



产品概述 (Summary)

产品名称 (Production Name)	CD315 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD315
别名 (Alternative Names)	PTGFRN; FPRP; EWI-F; CD9P-1; SMAP-6
基因 ID (Gene ID)	5738.0
蛋白 ID (SwissProt ID)	Q9P2B2.Purified recombinant fragment of human CD315 (AA: extra 46-221) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	98.6kDa

研究背景 (Background)

Inhibits the binding of prostaglandin F2-alpha (PGF2-alpha) to its specific FP receptor, by decreasing the receptor number


产品名称: CD315 Mouse Monoclonal Antibody
产品货号: AMM82296




rather than the affinity constant. Functional coupling with the prostaglandin F2-alpha receptor seems to occur (By similarity).

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD315 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD315 mouse mAb with DAB staining.

产品名称: CD315 Mouse Monoclonal Antibody
产品货号: AMM82296



注意事项 (Note)

For research use only .