

产品名称: CD369 Mouse Monoclonal Antibody  
产品货号: AMM82297



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | CD369 Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody       |
| 宿主 (Host)              | Mouse                           |
| 应用 (Application)       | ELISA,FC                        |
| 种属反应性 (Reactivity)     | Human                           |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide                                         |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                              |
|------------------------|----------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CD369                                                                                        |
| 别名 (Alternative Names) | CLEC7A; BGR; CANDF4; SCARE2; DECTIN1; CLECSF12                                               |
| 基因 ID (Gene ID)        | 64581.0                                                                                      |
| 蛋白 ID (SwissProt ID)   | Q9BXN2.Purified recombinant fragment of human CD369 (AA: extra 66-247) expressed in E. Coli. |

## 产品应用 (Application)

|                          |                                     |
|--------------------------|-------------------------------------|
| 稀释比 (Dilution Ratio)     | ELISA 1:5000-1:20000,FC 1:200-1:400 |
| 蛋白分子量 (Molecular Weight) | 27.6kDa                             |

## 研究背景 (Background)

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. The encoded


产品名称: CD369 Mouse Monoclonal Antibody  
产品货号: AMM82297




glycoprotein is a small type II membrane receptor with an extracellular C-type lectin-like domain fold and a cytoplasmic domain with an immunoreceptor tyrosine-based activation motif. It functions as a pattern-recognition receptor that recognizes a variety of beta-1,3-linked and beta-1,6-linked glucans from fungi and plants, and in this way plays a role in innate immune response. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.

## 研究领域 (Research Area)

## 图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD369 mouse mAb (green) and negative control (red).

## 注意事项 (Note)

For research use only .