

产品名称: CD307A Mouse Monoclonal Antibody
产品货号: AMM82300



产品概述 (Summary)

产品名称 (Production Name)	CD307A Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD307A
别名 (Alternative Names)	FCRL1; FCRH1; IFGP1; IRTA5
基因 ID (Gene ID)	115350.0
蛋白 ID (SwissProt ID)	Q96LA6.Purified recombinant fragment of human CD307A (AA: extra 17-202) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)

This gene encodes a member of the immunoglobulin receptor superfamily and is one of several Fc receptor-like


产品名称: CD307A Mouse Monoclonal Antibody
产品货号: AMM82300




glycoproteins clustered on the long arm of chromosome 1. The encoded protein contains three extracellular C2-like immunoglobulin domains, a transmembrane domain and a cytoplasmic domain with two immunoreceptor-tyrosine activation motifs. This protein may play a role in the regulation of cancer cell growth. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD307A mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .