

产品名称: RB1CC1 Mouse Monoclonal Antibody
产品货号: AMM82303



产品概述 (Summary)

产品名称 (Production Name)	RB1CC1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RB1CC1
别名 (Alternative Names)	CC1; ATG17; FIP200; PPP1R131
基因 ID (Gene ID)	9821.0
蛋白 ID (SwissProt ID)	Q8TDY2 .Purified recombinant fragment of human RB1CC1 (AA: 148-244) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	183kDa

研究背景 (Background)

The protein encoded by this gene interacts with signaling pathways to coordinately regulate cell growth, cell proliferation,


产品名称: RB1CC1 Mouse Monoclonal Antibody
产品货号: AMM82303



apoptosis, autophagy, and cell migration. This tumor suppressor also enhances retinoblastoma 1 gene expression in cancer cells. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)

注意事项 (Note)

For research use only .