

产品名称: PINK1 Mouse Monoclonal Antibody
产品货号: AMM82317



产品概述 (Summary)

产品名称 (Production Name)	PINK1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PINK1
别名 (Alternative Names)	BRPK; PARK6
基因 ID (Gene ID)	65018.0
蛋白 ID (SwissProt ID)	Q9BXM7.Purified recombinant fragment of human PINK1 (AA: 112-273) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	62.8kDa

研究背景 (Background)

This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from

产品名称: PINK1 Mouse Monoclonal Antibody
产品货号: AMM82317




stress-induced mitochondrial dysfunction. Mutations in this gene cause one form of autosomal recessive early-onset Parkinson disease.


研究领域 (Research Area)

Autophagy


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using PINK1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using PINK1 mouse mAb with DAB staining.

产品名称: PINK1 Mouse Monoclonal Antibody
产品货号: AMM82317



注意事项 (Note)

For research use only .