

产品名称: CD364 Mouse Monoclonal Antibody
产品货号: AMM82332



产品概述 (Summary)

产品名称 (Production Name)	CD364 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD364
别名 (Alternative Names)	PI16; PSPBP; CRISP9; MSMBBP
基因 ID (Gene ID)	221476.0
蛋白 ID (SwissProt ID)	Q6UXB8 .Purified recombinant fragment of human CD364 (AA: 79-210) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	49.5kDa

研究背景 (Background)

PI16 (Peptidase Inhibitor 16) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include


产品名称: CD364 Mouse Monoclonal Antibody
产品货号: AMM82332




peptidase inhibitor activity.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CD364 mouse mAb against THP-1 (1), PC-3 (2), LNCAP (3), Jurkat (4), and A431 (5) cell lysate.



Flow cytometric analysis of HL-60 cells using CD364 mouse mAb (green) and negative control (red).

产品名称: CD364 Mouse Monoclonal Antibody
产品货号: AMM82332



注意事项 (Note)

For research use only .