

产品名称: PMAIP1 Mouse Monoclonal Antibody
产品货号: AMM82353



产品概述 (Summary)

产品名称 (Production Name)	PMAIP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PMAIP1
别名 (Alternative Names)	APR; NOXA
基因 ID (Gene ID)	5366.0
蛋白 ID (SwissProt ID)	Q13794.Purified recombinant fragment of human PMAIP1 (AA: 1-54) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	6kDa

研究背景 (Background)

The protein encoded by this gene promotes activation of caspases and apoptosis, Promotes mitochondrial membrane

产品名称: PMAIP1 Mouse Monoclonal Antibody
产品货号: AMM82353




changes and efflux of apoptogenic proteins from the mitochondria. Contributes to p53/TP53-dependent apoptosis after radiation exposure. Promotes proteasomal degradation of MCL1. Competes with BAK1 for binding to MCL1 and can displace BAK1 from its binding site on MCL1 (By similarity). Competes with BIM/BCL2L11 for binding to MCL1 and can displace BIM/BCL2L11 from its binding site on MCL1.


研究领域 (Research Area)

Apoptosis

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using PMAIP1 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .