

产品名称: BRIP1 Mouse Monoclonal Antibody
产品货号: AMM82390



产品概述 (Summary)

产品名称 (Production Name)	BRIP1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BRIP1
别名 (Alternative Names)	OF; BACH1; FANCI
基因 ID (Gene ID)	83990.0
蛋白 ID (SwissProt ID)	Q9BX63.Purified recombinant fragment of human BRIP1 (AA: 904-986) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	140.8kDa

研究背景 (Background)

The protein encoded by this gene is a member of the RecQ DEAH helicase family and interacts with the BRCT repeats of


产品名称: BRIP1 Mouse Monoclonal Antibody
产品货号: AMM82390




breast cancer, type 1 (BRCA1). The bound complex is important in the normal double-strand break repair function of breast cancer, type 1 (BRCA1). This gene may be a target of germline cancer-inducing mutations.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of Hela cells using BRIP1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using BRIP1 mouse mAb with DAB staining.

产品名称: BRIP1 Mouse Monoclonal Antibody
产品货号: AMM82390



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using BRIP1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .