

产品名称: CD7 Mouse Monoclonal Antibody
产品货号: AMM82398



产品概述 (Summary)

产品名称 (Production Name)	CD7 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD7
别名 (Alternative Names)	GP40; TP41; Tp40; LEU-9
基因 ID (Gene ID)	924.0
蛋白 ID (SwissProt ID)	P09564.Purified recombinant fragment of human CD7 (AA: 26-180) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	25.4kDa

研究背景 (Background)

This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found


产品名称: CD7 Mouse Monoclonal Antibody
产品货号: AMM82398




on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of THP-1 cells using CD7 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using CD7 mouse mAb with DAB staining.

产品名称: CD7 Mouse Monoclonal Antibody
产品货号: AMM82398



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using CD7 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .