

产品名称: SOX11 Mouse Monoclonal Antibody
产品货号: AMM82420



产品概述 (Summary)

产品名称 (Production Name)	SOX11 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SOX11
别名 (Alternative Names)	CSS9; MRD27
基因 ID (Gene ID)	6664.0
蛋白 ID (SwissProt ID)	P35716.Purified recombinant fragment of human SOX11 (AA: 1-250) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	46.7kDa

研究背景 (Background)

This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the


产品名称: SOX11 Mouse Monoclonal Antibody
产品货号: AMM82420




regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The protein may function in the developing nervous system and play a role in tumorigenesis.

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Hela cells using SOX11 mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .