

产品名称: Vimentin Mouse Monoclonal Antibody
产品货号: AMM82438



产品概述 (Summary)

产品名称 (Production Name)	Vimentin Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Vimentin
别名 (Alternative Names)	FLJ36605; VIM
基因 ID (Gene ID)	7431.0
蛋白 ID (SwissProt ID)	P08670.Purified recombinant fragment of human VIM (AA: 2-466) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	53.7kDa

研究背景 (Background)

This gene encodes a type III intermediate filament protein. Intermediate filaments, along with microtubules and actin


产品名称: Vimentin Mouse Monoclonal Antibody
产品货号: AMM82438




microfilaments, make up the cytoskeleton. The encoded protein is responsible for maintaining cell shape and integrity of the cytoplasm, and stabilizing cytoskeletal interactions. This protein is involved in neuritogenesis and cholesterol transport and functions as an organizer of a number of other critical proteins involved in cell attachment, migration, and signaling. Bacterial and viral pathogens have been shown to attach to this protein on the host cell surface. Mutations in this gene are associated with congenital cataracts in human patients.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using VIM mouse mAb against Jurkat (1), K562 (2), SK-N-SH (3), SH-SY5Y (4), HeLa (5), NIH/3T3 (6), C6 (7), and RAW264.7 (8) cell lysate.


产品名称: Vimentin Mouse Monoclonal Antibody
产品货号: AMM82438



Flow cytometric analysis of Hela cells using VIM mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using VIM mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using VIM mouse mAb with DAB staining.

注意事项 (Note)

For research use only .