

产品名称: MLH1 Mouse Monoclonal Antibody
产品货号: AMM82444



产品概述 (Summary)

产品名称 (Production Name)	MLH1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLH1
别名 (Alternative Names)	FCC2; COCA2; HNPCC; hMLH1; HNPCC2
基因 ID (Gene ID)	4292.0
蛋白 ID (SwissProt ID)	P40692.Purified recombinant fragment of human MLH1 (AA:381-483) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	84.6kDa

研究背景 (Background)

The protein encoded by this gene can heterodimerize with mismatch repair endonuclease PMS2 to form MutL alpha, part


产品名称: MLH1 Mouse Monoclonal Antibody
产品货号: AMM82444




of the DNA mismatch repair system. When MutL alpha is bound by MutS beta and some accessory proteins, the PMS2 subunit of MutL alpha introduces a single-strand break near DNA mismatches, providing an entry point for exonuclease degradation. The encoded protein is also involved in DNA damage signaling and can heterodimerize with DNA mismatch repair protein MLH3 to form MutL gamma, which is involved in meiosis. This gene was identified as a locus frequently mutated in hereditary nonpolyposis colon cancer (HNPCC). [provided by RefSeq, Aug 2017]

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using MLH1 mouse mAb against Hela (1), Jurkat (2), A431 (3), HepG2 (4), and MCF-7 (5) cell lysate.


产品名称: MLH1 Mouse Monoclonal Antibody
产品货号: AMM82444



Flow cytometric analysis of Hela cells using MLH1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human Colon cancer tissues using MLH1 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human Ovarian cancer tissues using MLH1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .