

产品名称: WT1 Mouse Monoclonal Antibody
产品货号: AMM82467



产品概述 (Summary)

产品名称 (Production Name)	WT1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,ELISA,FC
种属反应性 (Reactivity)	Human, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WT1
别名 (Alternative Names)	GUD; AWT1; WAGR; WT33; NPHS4; WIT-2
基因 ID (Gene ID)	7490.0
蛋白 ID (SwissProt ID)	P19544.Purified recombinant fragment of human WT1 (AA: 1-181) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	49.2kDa

研究背景 (Background)


产品名称: WT1 Mouse Monoclonal Antibody
产品货号: AMM82467




This gene encodes a transcription factor that contains four zinc-finger motifs at the C-terminus and a proline/glutamine-rich DNA-binding domain at the N-terminus. It has an essential role in the normal development of the urogenital system, and it is mutated in a small subset of patients with Wilms tumor. This gene exhibits complex tissue-specific and polymorphic imprinting pattern, with biallelic, and monoallelic expression from the maternal and paternal alleles in different tissues. Multiple transcript variants have been described. In several variants, there is evidence for the use of a non-AUG (CUG) translation initiation codon upstream of, and in-frame with the first AUG. Authors of PMID:7926762 also provide evidence that WT1 mRNA undergoes RNA editing in human and rat, and that this process is tissue-restricted and developmentally regulated.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using WT1 mouse mAb against HEK293 (1), COS7 (2), and PC-3 (3) cell lysate.


产品名称: WT1 Mouse Monoclonal Antibody
产品货号: AMM82467



Immunofluorescence analysis of HeLa cells using WT1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)




Flow cytometric analysis of HeLa cells using WT1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using WT1 mouse mAb with DAB staining.

产品名称: WT1 Mouse Monoclonal Antibody
产品货号: AMM82467



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using WT1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .