

产品名称: CD267 Mouse Monoclonal Antibody
产品货号: AMM82486



产品概述 (Summary)

产品名称 (Production Name)	CD267 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD267
别名 (Alternative Names)	CVID; RYZN; TACI; CD267; CVID2; IGAD2; TNFRSF14B
基因 ID (Gene ID)	23495.0
蛋白 ID (SwissProt ID)	O14836.Purified recombinant fragment of human CD267 (AA: extra(1-165)) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	31.8kDa

研究背景 (Background)

The protein encoded by this gene is a lymphocyte-specific member of the tumor necrosis factor (TNF) receptor superfamily.


产品名称: CD267 Mouse Monoclonal Antibody
产品货号: AMM82486




It interacts with calcium-modulator and cyclophilin ligand (CAML). The protein induces activation of the transcription factors NFAT, AP1, and NF-kappa-B and plays a crucial role in humoral immunity by interacting with a TNF ligand. This gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of Jurkat cells using CD267 mouse mAb (green) and negative control (red).




Flow cytometric analysis of THP-1 cells using CD267 mouse mAb (green) and negative control (red).

产品名称: CD267 Mouse Monoclonal Antibody
产品货号: AMM82486



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD267 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human lymphoid tissues using CD267 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .