

产品名称: CHGA Mouse Monoclonal Antibody
产品货号: AMM82506



产品概述 (Summary)

产品名称 (Production Name)	CHGA Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHGA
别名 (Alternative Names)	CGA
基因 ID (Gene ID)	1113.0
蛋白 ID (SwissProt ID)	P10645.Purified recombinant fragment of human CHGA (AA: 278-457) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	50kDa

研究背景 (Background)

The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory


产品名称: CHGA Mouse Monoclonal Antibody
产品货号: AMM82506




proteins. It is found in secretory vesicles of neurons and endocrine cells. This gene product is a precursor to three biologically active peptides; vasostatin, pancreastatin, and parastatin. These peptides act as autocrine or paracrine negative modulators of the neuroendocrine system. Two other peptides, catestatin and chromofungin, have antimicrobial activity and antifungal activity, respectively. Two transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using CHGA mouse mAb against SH-SY5Y (1) cell lysate.




产品名称: CHGA Mouse Monoclonal Antibody
产品货号: AMM82506



Flow cytometric analysis of C6 cells using CHGA mouse mAb (green) and negative control (red).



Flow cytometric analysis of Hela cells using CHGA mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human pancreatic tissues using CHGA mouse mAb with DAB staining.

注意事项 (Note)

For research use only .