

产品名称: GLUL Mouse Monoclonal Antibody
产品货号: AMM82510



产品概述 (Summary)

产品名称 (Production Name)	GLUL Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GLUL
别名 (Alternative Names)	GS; GLNS; PIG43; PIG59
基因 ID (Gene ID)	2752.0
蛋白 ID (SwissProt ID)	P15104.Purified recombinant fragment of human GLUL (AA: 2-121) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	42kDa

研究背景 (Background)

The protein encoded by this gene belongs to the glutamine synthetase family. It catalyzes the synthesis of glutamine from


产品名称: GLUL Mouse Monoclonal Antibody
产品货号: AMM82510




glutamate and ammonia in an ATP-dependent reaction. This protein plays a role in ammonia and glutamate detoxification, acid-base homeostasis, cell signaling, and cell proliferation. Glutamine is an abundant amino acid, and is important to the biosynthesis of several amino acids, pyrimidines, and purines. Mutations in this gene are associated with congenital glutamine deficiency, and overexpression of this gene was observed in some primary liver cancer samples. There are six pseudogenes of this gene found on chromosomes 2, 5, 9, 11, and 12. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using GLUL mouse mAb against Jurkat (1), and mouse liver (2) cell lysate.


产品名称: GLUL Mouse Monoclonal Antibody
产品货号: AMM82510



Flow cytometric analysis of Hela cells using GLUL mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebellar tissues using GLUL mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GLUL mouse mAb with DAB staining.

注意事项 (Note)

For research use only .