

产品名称: SALL4 Mouse Monoclonal Antibody
产品货号: AMM82514



产品概述 (Summary)

产品名称 (Production Name)	SALL4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA,FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SALL4
别名 (Alternative Names)	DRRS; HSAL4; ZNF797; dJ1112F19.1
基因 ID (Gene ID)	57167.0
蛋白 ID (SwissProt ID)	Q9UJQ4.Purified recombinant fragment of human SALL4 (AA: 96-359) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	112kDa

研究背景 (Background)

The protein encoded by this gene may be a zinc finger transcription factor. Defects in this gene are a cause of Duane-radial


产品名称: SALL4 Mouse Monoclonal Antibody
产品货号: AMM82514




ray syndrome (DRRS).

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using SALL4 mouse mAb against NIH/3T3 (1) C2C12 (2) cell lysate.



Flow cytometric analysis of Hela cells using SALL4 mouse mAb (green) and negative control (red).

产品名称: SALL4 Mouse Monoclonal Antibody
产品货号: AMM82514



注意事项 (Note)

For research use only .