

产品名称: 5-hydroxyuridine Mouse Monoclonal Antibody
产品货号: AMM82518



产品概述 (Summary)

产品名称 (Production Name)	5-hydroxyuridine Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	5-hydroxyuridine
别名 (Alternative Names)	5-hydroxyuridine
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	/Synthesized peptide of 5-hydroxyuridine-KLH

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa


研究背景 (Background)

产品名称: 5-hydroxyuridine Mouse Monoclonal Antibody
产品货号: AMM82518




研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using 5-hydroxyuridine mouse mAb against mouse serum (1) and rat serum (2) lysate.

注意事项 (Note)

For research use only .