

产品名称: FLI1 Mouse Monoclonal Antibody
产品货号: AMM82520



产品概述 (Summary)

产品名称 (Production Name)	FLI1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FLI1
别名 (Alternative Names)	EWSR2; SIC-1; BDPLT21
基因 ID (Gene ID)	2313.0
蛋白 ID (SwissProt ID)	Q01543.Purified recombinant fragment of human FLI1 (AA: 303-452) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	50.9kDa

研究背景 (Background)

This gene encodes a transcription factor containing an ETS DNA-binding domain. The gene can undergo a t(11;22)


产品名称: FLI1 Mouse Monoclonal Antibody
产品货号: AMM82520




(q24;q12) translocation with the Ewing sarcoma gene on chromosome 22, which results in a fusion gene that is present in the majority of Ewing sarcoma cases. An acute lymphoblastic leukemia-associated t(4;11)(q21;q23) translocation involving this gene has also been identified. Alternative splicing results in multiple transcript variants.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of Hela cells using FLI1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human breast cancer tissues using FLI1 mouse mAb with DAB staining.

产品名称: FLI1 Mouse Monoclonal Antibody
产品货号: AMM82520



Immunohistochemical analysis of paraffin-embedded human kidney tissues using FLI1 mouse mAb with DAB staining.

注意事项 (Note)

For research use only .