

产品名称: CD15 Mouse Monoclonal Antibody
产品货号: AMM82534



产品概述 (Summary)

产品名称 (Production Name)	CD15 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD15
别名 (Alternative Names)	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
基因 ID (Gene ID)	2526.0
蛋白 ID (SwissProt ID)	P22083.Purified recombinant fragment of human CD15 (AA: 1-147) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	59kDa

研究背景 (Background)

The product of this gene transfers fucose to N-acetyllactosamine polysaccharides to generate fucosylated carbohydrate


产品名称: CD15 Mouse Monoclonal Antibody
产品货号: AMM82534




structures. It catalyzes the synthesis of the non-sialylated antigen, Lewis x (CD15). [provided by RefSeq, Jan 2009]

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using CD15 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using CD15 mouse mAb (green) and negative control (red).

产品名称: CD15 Mouse Monoclonal Antibody
产品货号: AMM82534



注意事项 (Note)

For research use only .