

产品名称: CD21 Mouse Monoclonal Antibody
产品货号: AMM82540



产品概述 (Summary)

产品名称 (Production Name)	CD21 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD21
别名 (Alternative Names)	CR2; CR; C3DR; CVID7; SLEB9
基因 ID (Gene ID)	1380.0
蛋白 ID (SwissProt ID)	P20023.Purified recombinant fragment of human CD21 (AA: extra 740-964) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	112.9kDa

研究背景 (Background)

This gene encodes a membrane protein, which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T


产品名称: CD21 Mouse Monoclonal Antibody
产品货号: AMM82540




lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

研究领域 (Research Area)


图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD21 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using CD21 mouse mAb (green) and negative control (red).

产品名称: CD21 Mouse Monoclonal Antibody
产品货号: AMM82540



注意事项 (Note)

For research use only .