

产品名称: C-MYC Mouse Monoclonal Antibody
产品货号: AMM82562



产品概述 (Summary)

产品名称 (Production Name)	C-MYC Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC,ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C-MYC
别名 (Alternative Names)	MRTL; MYCC; c-Myc; bHLHe39
基因 ID (Gene ID)	4609.0
蛋白 ID (SwissProt ID)	P01106.Purified recombinant fragment of human C-MYC (AA: (290-439)) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	32.1/57/48.8kDa

研究背景 (Background)

This gene is a proto-oncogene and encodes a nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis

产品名称: C-MYC Mouse Monoclonal Antibody
产品货号: AMM82562




and cellular transformation. The encoded protein forms a heterodimer with the related transcription factor MAX. This complex binds to the E box DNA consensus sequence and regulates the transcription of specific target genes. Amplification of this gene is frequently observed in numerous human cancers. Translocations involving this gene are associated with Burkitt lymphoma and multiple myeloma in human patients. There is evidence to show that translation initiates both from an upstream, in-frame non-AUG (CUG) and a downstream AUG start site, resulting in the production of two isoforms with distinct N-termini. [provided by RefSeq, Aug 2017]


研究领域 (Research Area)

PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of HeLa cells using C-MYC mouse mAb (green) and negative control (red).


产品名称: C-MYC Mouse Monoclonal Antibody
产品货号: AMM82562



Flow cytometric analysis of Hepg2 cells using C-MYC mouse mAb (green) and negative control (red).




Flow cytometric analysis of Lovo cells using C-MYC mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using C-MYC mouse mAb with DAB staining.

产品名称: C-MYC Mouse Monoclonal Antibody
产品货号: AMM82562



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using C-MYC mouse mAb with DAB staining.

注意事项 (Note)

For research use only .