

产品名称: CD62E Mouse Monoclonal Antibody
产品货号: AMM82600



产品概述 (Summary)

产品名称 (Production Name)	CD62E Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	ELISA,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD62E
别名 (Alternative Names)	ELAM; ESEL; SELE; ELAM1; LECAM2
基因 ID (Gene ID)	6401.0
蛋白 ID (SwissProt ID)	P16581.Purified recombinant fragment of human CD62E (AA: extra(22-162)) expressed in E. Coli.

产品应用 (Application)

稀释比 (Dilution Ratio)	ELISA 1:5000-1:20000,FC 1:200-1:400
蛋白分子量 (Molecular Weight)	66.7kDa

研究背景 (Background)

The protein encoded by this gene is found in cytokine-stimulated endothelial cells and is thought to be responsible for the


产品名称: CD62E Mouse Monoclonal Antibody
产品货号: AMM82600




accumulation of blood leukocytes at sites of inflammation by mediating the adhesion of cells to the vascular lining. It exhibits structural features such as the presence of lectin- and EGF-like domains followed by short consensus repeat (SCR) domains that contain 6 conserved cysteine residues. These proteins are part of the selectin family of cell adhesion molecules. Adhesion molecules participate in the interaction between leukocytes and the endothelium and appear to be involved in the pathogenesis of atherosclerosis. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using CD62E mouse mAb (green) and negative control (red).

注意事项 (Note)

For research use only .